WSB

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Contents

1 Introduction	13
2 Safety Instructions	13
2.1 Description of symbols	13
2.2 Instructions	14
3 Accessoires	15
4 Product description	15
4.1 Key features	15
4.2 General survey	16
4.3 Operating surface	16
4.4 Term description	16
5 Positioning the device and commissioning.....	17
6 Job description.....	18
6.1 Menu structure.....	18
6.2 Set turn-off delay	18
6.3 Set turn-on delay.....	18
6.4 Set temperature.....	19
6.5 Set shaking speed	19
6.6 Justification.....	19
6.7 Key lock.....	19
6.8 Storage function	19
7 Maintenance and cleaning.....	20
8 Troubleshooting Guide	20
9 Technical data.....	21
10 Certificate.....	22

1 Introduction

Thank you for purchasing our **WSB Shaking Water Bath**. You did a good choice buying a **Wisd** product.

Before use, please read this manual and use it properly.

To cover the user's health and avoid the product's destruction you have to follow the safety instructions of the second chapter.

Please keep this operating manual in a save place for future reference and hand it out if you consign the device to a third person.

Witeg reserves to modify the capsule or make amendments for the improvement of quality and efficiency without prior notice.

2 Safety Instructions

2.1 Description of symbols

**WARNING**

Information for avoiding injuries or fatal accidents.

**ATTENTION**

Safety instruction to avoid a damage of the product.

**ATTENTION**

Information about technical parameters.

**Environment**

Information about the operating environment.

**WARNING of explosions**

Safety instruction to avoid explosions.

**ATTENTION to high temperatures**

Instructions to avoid injuries caused by high temperatures.

**Risk of injury**

Avoid injuries of your hands by following these instructions.

2.2 Instructions

For a safety and faultless using of the device it is indispensable that you read the following instructions purposely and follow them by the use.

- Connect the device only to the suitable power supply. Please have a look to the identification plate.
- Use the device only in electric circuits with ground fault circuit interrupter.
- Check the power supply cord before you stick it into the electrical outlet.
- Do not use power bars to extend the power supply cord.
- Only use extension cables with sufficient cable cross section.
- Do not touch the device with wet hands during the usage because this might cause an electrical shock.
- If unexpected sound, smell or smog is generated by the device, pull out the main plug and contact the manufacturer or your supplier immediately.
- Disconnect the power supply plug before cleaning or maintenance.
- If you use radiation or contaminated samples the responsibility is only up to you.

- It is important that the device stands up on a horizontal, stabile and firm subsurface to avoid vibrations.
- Make sure that the device cannot slip away.
- To avoid heat congestion or fires make sure that there is enough space for air circulations.
- Avoid direct solar radiation.
- Do not use the device near highly combustible materials.
- Do not set the device outdoors.
- It is superiorly to set the device in a good illuminated place to avoid accidents engendered by handling errors.
- To set the device in wet or dusty places could cause overheating, short circuits or fires.

- Do not use this device in area with potentially explosive atmosphere.

- Reparations and amendments should only made by **witeg** authorized people. Contraventions exclude the warranty.
- Make sure, that the safety switch is adjusted up to 30°C over max. temperature. Otherwise the heating maybe doesn't work.

- Do not put heat-sensitive items at the device.
- Warn all attendees if you leave the device unattended.
- Pay attention when you lift up the lid, hot condensation water may drop down or hot steam could come out.
- Do not clean the scorching device.

- Use this device only if you read the whole safety instructions.
- This device should only use by persons who are familiarized with safety instructions made for laboratories.
- Pay attention that no liquid reaches the controller or inward of the capsule.
- You have to keep the device clean if you want to get a longer life of it.
- Wear corresponding protective work clothing during the usage.
- Pull out the power cable if you didn't want to use the device for a longer time.

- Do not put your hands between shaking rack and bath during the shaking process.
- Pull out the power supply cord before you clean the device.

3 Accessoires

Quantity	DH.WSB03018	DH.WSB03030	DH.WSB03045
1	WSB-18	WSB-30	WSB-45
1	Power cord		
1	Operating manual		

4 Product description

4.1 Key features

The **Wisdom WSB Shaking Water Bath** is designed for chemical, biological, medical and pharmaceutical applications. The device was developed for simplest usage and a very long life.

Key features are as follows:


- Simplest operating by MAIN button and Jog-Dial knob.
- Shaking rack for flasks and tubes up to 2000 ml.
- Fast commissioning.
- Long life by simply and solid construction.
- Highest temperature accuracy by digital fuzzy controller.
- Temperature, speed and delay time showed on illuminated LCD.
- Storage functions for set time and temperature values.
- Turn-off and -on delay timer up to 99 h 59 min.

- Signal sound after end of turn-off delay.
- Easy to justify.
- Key lock against inadvertent changing of set speed and timer values.
- Over temperature protection by safety switch.

4.2 General survey

4.3 Operating surface

4.4 Term description

SV: set value (upper line, in key lock menu no SV)

PV: present value, measured by the device (lower line, in key lock menu also upper line)

Main switch

- On / off switch.

Timer lamp

- Is blinking in menu items for timer modulation
- Shines continuous in the key lock menu item.

Heating lamp

- Shines if the heater is activated.

Shaking lamp

- Shines if shaking is activated.

Operating lamp

- Shines if heater, timer or shaker is activated.

Jog-Dial knob

- Changes values by turning and stores, confirms, stops and starts functions by pushing.

Safety switch

- The safety switch turns the heating off if your set maximal temperature is overshoot.
- Set the maximal temperature roughly 30 °C over the desired maximal temperature, otherwise the device maybe not reaches your set value.

MAIN button

- For changing between the different menu items.

SUB button

- To reach the sub menu items.


5 Positioning the device and commissioning

1. Follow especially the safety instructions marked with following symbol:
2. Remove the package and perhaps soiling carefully.
3. If there are any unexpected transport damages, please inform your forwarder and your customer service.
4. It is important that the device stands up on a horizontal, stabile and firm subsurface to avoid vibrations.
5. Set the protection switch 30 °C over desired bath temperature otherwise the heating coil does not work.
6. Fill enough distilled water in.
7. Make sure that you follow all safety instructions.
8. Before you connect the device with the electrical outlet, please let some time laps until the device is acclimatized.
9. Connect the device with the power outlet.
10. Switch the device on and start the process (look at job description).
11. Make sure that the device cannot slip away.

6 Job description

6.1 Menu structure

6.2 Set turn-off delay

The turn-off delay timer serves as automatically stop function. It turns off the heating and shaking after the timeframe you set is elapsed.

- Use MAIN button to reach next menu without setting a delay time.
- The delay time is adjustable from one minute till over 99 hours (00:01 - 99:59) in one minute steps.
- Turn the Jog-Dial knob to set your delay time. Push the Jog-Dial if you reached the desired value to confirm it. With the confirmation the timeframe starts elapsing.
- If you want to change the value after confirmation, you have to push the Jog-Dial once again. The PV initiates blinking and you can modify the value again.
- To confirm the value and start the delay time push the Jog-Dial yet again.

6.3 Set turn-on delay

The turn-on delay timer serves as automatically start function. It turns on the heating after the timeframe you set is elapsed.

To reach this sub menu item you have to use the SUB button in the menu “set turn-off delay”

- The delay time is adjustable from one minute till over 99 hours (00:01 - 99:59) in one minute steps.
- Turn the Jog-Dial knob to set your delay time. Push the Jog-Dial if you reached the desired value to confirm it. With the confirmation the timeframe starts elapsing.

- If you want to change the value after confirmation, you have to push the Jog-Dial once again. The PV initiates blinking and you can modify the value again.
- To confirm the value and start the delay time push the Jog-Dial yet again.

6.4 Set temperature

- When you push the MAIN button in the menu item “set turn-off delay” you reach the menu item “set temperature”. You can recognize this menu at “°C” in SV and PV line.
- Push the MAIN button once again if you want to work without the heating function. Then you reach the next menu without a collimated temperature value.
- Choose your desired temperature by turning the Jog-Dial and confirm the set value by pushing the Jog-Dial.
- For changing the confirmed value again, you first have to push the Jog-Dial knob again. “HEATING” starts blinking at the LCD and the temperature is adjustable again.
- When the heating lamp is shining, the heater is activated.

6.5 Set shaking speed

- When you push the MAIN button in the menu item “set temperature” you reach the menu item “set shaking speed”.
- Choose your desired speed by turning the Jog-Dial and confirm the set value by pushing the Jog-Dial.
- For changing the confirmed value again, you first have to push the Jog-Dial knob again. SV starts blinking at the LCD and the speed is adjustable again.

6.6 Justification

If the temperature measured by you differs from the temperature the device is showing (PV) you have the possibility to compensate the divergence by the Offset. So you can balance differences of about ± 10 °C.

The Offset is calculated as follows:

Temperature measured by you – showed temperature (PV) = Offset

- To reach this menu you have to use the SUB button in the menu item “set temperature”.
- Set the calculated value by turning the Jog-Dial. The Offset will be added to the PV immediately.
- Save the offset by holding the Jog-Dial pushed until the values in PV and SV line starts blinking.
- Push the MAIN button to leave the justification menu.

6.7 Key lock

- To avoid accidental changing of the set values, push the MAIN button in the menu item “set temperature”.
- In this menu item can you see the actually temperature in the SV-line and the remaining time of turn-off delay in the PV-line. The function of the Jog-Dial knob is deactivated.
- Press the MAIN button once more to leave the key lock. Than you are in the menu item “set turn-off delay”.

6.8 Storage function

The **WSB Shaking Water Bath** includes a storage function for temperature, speed and timer values. That is practical if you want to work with the same parameters often. Then you did not

need to set them in every usage again.

- For storing values you only need to hold the Jog-Dial knob pushed for two seconds. Successful storing is displayed by blinking of the SV for three times.
- Any stored parameters are erasable. Turn the value with the Jug-Shuttle till "0" and then hold the knob pressed for two seconds.
- If you want to overwrite a value it is not necessary to erase the old one. Set the desired value and use the Jog-Dial knob for storing.

7 Maintenance and cleaning

- Please read the whole safety instructions before you start to maintain or clean the device.
- Pull the power plug out of the electrical outlet.
- Only use care cleaning agents and soft cleaning rags.
- Organic lotions, strong chemicals and rough cleaning rags can damage the device.
- If you not want to use the device for a long time, pull of the power plug and keep it dry in the package.
- Only use original spare parts.

8 Troubleshooting Guide

Situation	Confirmations and Solutions
Power Failure	<ul style="list-style-type: none"> ✓ Check the supplied voltage is present. ✓ Check the power connection and the power cord for proper connection and integrity. ✓ Check if the MAIN Power Switch is "ON". ✓ Fuse is blown.
Unit Not Heating	<ul style="list-style-type: none"> ✓ Check if the MAIN Power Switch is "ON". ✓ Check if the set temperature (SV) confirmed with the Jog-Dial knob (by turn the Jog-Dial knob cannot change the set temperature – see point 6.4) ✓ Check the Safety switch is set higher (min. 10°C-20°C) than the set temperature (SV).
Heating, but unit not reaching temperature	<ul style="list-style-type: none"> ✓ Check the Safety switch is set higher (min. 10°C-20°C) than the set temperature (SV). ✓ Problem in the sensor or main controller.
Unit Not Shaking	<ul style="list-style-type: none"> ✓ Check if the MAIN Power Switch is "ON". ✓ Check if the set RPM confirmed with the Jog-Dial knob (see point 6.5).
ERROR 1	<ul style="list-style-type: none"> ✓ Temperature Sensor Error.

If other problems arise, or one of the above not solves problems, please contact your official agent or the manufacturer.

Reparations and amendments should only made by **witeg** authorized people. Contraventions exclude the warranty.

9 Technical data

		WSB-18	WSB-30	WSB-45
Dimensions (B x T x H) (cm)	Bath	300 x 300 x 200	500 x 300 x 200	600 x 380 x 200
	Capsule	563 x 353 x 340	763 x 353 x 340	863 x 433 x 340
Capacity		18 l	30 l	45 l
Power consumption		1,4 kW	1,8 kW	2,4 kW
Capacity shaking rack		8 x 250 ml	15 x 250 ml	28 x 250 ml
Temp.	Range	Ambient temperature + 5 °C - 100 °C		
	Accuracy	± 0,1 °C		
	Sensor	PT100		
Speed		20 - 250 / min		
Controller		Digital Fuzzy-Controller with Jog-Dial knob (Turn +Push)		
Display		Illuminated LCD		
Timer		99 h 59 min (for Turn-off and turn-on delay timer)		
Material	Bath	Stainless steel		
	Capsule	Powder coated steel		
Other		<ul style="list-style-type: none"> - Storage function (Temp., Timer, Speed) - Key lock (Jog-Dial deactivated) - Easy to justify - Signal sound for timer 		
Power supply		AC 230 V, 1N~, 50/60 Hz		

10 Certificate

Certificate

for

The quality and all features were checked by the manufacturer before the shipment.

We grant from date of purchase

two years guaranty.

This certificate excludes damages by natural disasters or incorrect usages by the costumer.

Please look on your account and complete following table:

Article	WSB Shaking Water Bath
Typ	
Serialno.	
Date	

witeg Labortechnik GmbH

Am Bildacker 16
D-97877 Wertheim
TEL: +49-9342-9301-0
FAX: +49-9342-9301-77
Email: info@witeg.de
www.witeg.de