

# MTV-100 Multi-Tube Vortex



Programmable pulsing mode allows users to adjust the pulse-on and pulse-off times between 1~59 seconds in 1 second intervals. This feature enhances the vortex action by creating a more vigorous mix

## **Pulsing Mode:**

Programmable pulsing mode allows users to adjust the pulse-on and pulse-off times between 1~59 seconds in 1 second intervals. This feature enhances the vortex action by creating a more vigorous mix

## **Microprocessor Control:**

The variable speed microprocessor control provides consistent uniform mixing action

## **LCD Display:**

Independent LCD displays for speed and time allow operator to view both settings

## Applications:

Suspensions, high throughput testing in clinical, environmental, and chemistry labs

Model	MTV-100
Mixing Orbit	3.6mm
Speed Range	500~2500rpm
Timer	1min~99h59min
Speed Accuracy	10rpm
Maximum Weight Capacity	4.5kg
Tray Dimension(L X W)	310 x 180mm
Overall Dimension(Lx W x H)	426 x 250 x 480mm
Weight	21.5kg
Power	60W

## Accessories

Code	Description		
AS-02040-00	MTV-100 AC100~240V, 1.5A		
Code	Type	Tube Capacity	Dimension (mm)
AS-02041-01	MTV-B01	50 x $\Phi$ 10mm tube racks	245x132x45
AS-02041-02	MTV-B02	50 x $\Phi$ 12mm tube racks	245x132x45
AS-02041-03	MTV-B03	50 x $\Phi$ 13mm tube racks	245x132x45
AS-02041-04	MTV-B04	50 x $\Phi$ 16mm/15ml tube racks	245x132x45
AS-02041-05	MTV-B05	15 x $\Phi$ 25mm tube racks	245x132x45
AS-02041-06	MTV-B06	15 x $\Phi$ 29mm/50ml tube racks	245x132x45
AS-02041-00	MTV-B00	Tray pad set (upper & lower)	305×178.5×25