



YASUDA
— TESTING THE LIMITS —

No.546 Automatic Hoover Muller【YASUDA SEIKI】To grind・scatter the colorant

No.546 AUTOMATIC HOOVER MULLER



JIS-K5101

- This machine is used to grind and scatter colorant according to the Hoover Muller's method.

- After grinding and scattering the colorant, the operator is to use the GRIND METER to check the grind level of the test specimen.

Specification

Specifications are subject to change without notice.	
Glass Plate	ϕ 180 mm, Upper and Lower 1 pc each
Plate Rotation Speed	100 rpm
Weight Load	50 lbf, 100 lbf, 150 lbf [50 lbf = 22.68 kgf]
Counter	6 Digits Preset Counter, 6 Digits Running Counter
Power Source	AC 100 V, 1-Phase, 15 A, 50/60 Hz
Dimensions/ Weight (Approx.)	W500 × D500 × H550 mm, 55 kg