

D LAB orbital and linear shakers are ideal for almost any vessel from microcentrifuge tubes through petri dishes and microtitre plates to conical flasks.

SK-O330-Pro **SK-O180-Pro** LCD Digital Orbital Shaker

Features

- Max. load capacity is 7.5kg(O330) and 2.5kg(O180) with platform
- Double LCDs for independent display of time and speed
- Used with a wide variety of accessories
- Maintenance- free brushless DC motor
- Over speed detection and protection
- Remote function provides PC control and data transmission



Timer



Wide range of accessories

SK-O180-S LED Digital Orbital Shaker

Features

- Max. load capacity is 3kg with platform
- Double LEDs for independent display of time and speed
- Speed range of 40-200rpm
- Wide range of platforms for use with a variety of vessels






Specifications

	SK-O330-Pro	SK-O180-Pro	SK-O180-S
Orbital diameter	10mm	4mm	20mm
Max. load capacity (with platform)	7.5kg	2.5kg	3kg
Motor type	Brushless DC motor	Brushless DC motor	DC motor
Motor input power	28W	28W	16W
Motor output power	15W	15W	10W
Speed range	100-500rpm	100-800rpm	40-200rpm
Speed display	LCD	LCD	LED
Timer display	LCD	LCD	LED
Time settings range	1min-19h59min	1min-19h59min	1min-19h59min
Voltage	100-240V,50/60Hz	100-240V,50/60Hz	100-240V,50/60Hz
Power	30W	30W	20W
Weight	13.5kg	8.1kg	3.1kg
Permissible ambient temperature and humidity	5-40°C, 80%RH	5-40°C, 80%RH	5-40°C,80%RH
Dimension[W×D×H]	370×420×100mm	300×340×100mm	270×330×110mm
Protection class	IP21	IP21	IP21
Data connector	RS232	RS232	-





SK330.5



SK180.4



SK180.5



SK180.6



SK330.1



SK330.3



SK180.1



SK180.3



Cat. No.

Description

SK-O330-Pro
SK-L330-Pro
SK-R330-Pro
SK-D3309-Pro

SK-O180-Pro/S
SK-L180-Pro/S
SK-R1807-S
SK-D1807-S

18900155

Dish platform with non-slip mat and fixed string (8pcs)
Dimension[WxD]: 30x28.5cm



18202835

Silicone mattress for 16x 5~15ml tubes, use with SK330.5



18202846

Silicone mattress for 6x 50ml tubes, use with SK330.5



18900185

Dish attachment with non-slip mat
Dimension[WxD]: 24x21.5cm



18900200

Large size dish attachment with non-slip mat
Dimension[WxD]: 30x23cm



18900326

Dish attachment with non-slip mat
Dimension[WxD]: 28x28cm



18202852

Silicone mattress for 10x 5~10ml tubes



18202853

Silicone mattress for 4x 50ml tubes



18900027

Universal attachment with 4 bars
Size: 32x32cm
Capacity for flask: 25mlx21pcs,
50mlx15pcs, 100mlx12pcs,
250mlx6pcs, 500x4pcs



18900040

Lengthways roller attachment
Size: 32x32cm
Capacity for flask: 25mlx21pcs,
50mlx15pcs, 100mlx12pcs,
250mlx6pcs, 500x4pcs



18900025

Universal attachment with 3 bars
Size: 25x24.5cm
Capacity for flask: 25mlx10pcs,
50mlx8pcs, 100mlx6pcs,
250mlx4pcs, 500x2pcs



18900038

Lengthways roller attachment
Size: 25x24.5cm
Capacity for flask: 25mlx10pcs,
50mlx8pcs, 100mlx6pcs,
250mlx4pcs, 500x2pcs



18900036

Bar for SK330.1 universal attachment



18900037


Bar for SK180.1 universal attachment




18900280

Connect sticks for double dishes shaker





SK330.2



SK180.2



Cat. No.	Description	SK-O330-Pro SK-L330-Pro SK-R330-Pro SK-D3309-Pro	SK-O180-Pro/S SK-L180-Pro/S SK-R1807-S SK-D1807-S
18900028	Fixing clip attachment Size: 32 x 32cm	■	
18900026	Fixing clip attachment Size: 25 x 24cm		■
18900029	Fixing clip for flask volume 25mL, used with SK330.2[56pcs] SK180.2[36pcs]	■	■
18900030	Fixing clip for flask volume 50mL, used with SK330.2[28pcs] SK180.2[18pcs]	■	■
18900031	Fixing clip for flask volume 100mL, used with SK330.2[28pcs] SK180.2 [18pcs]	■	■
18900032	Fixing clip for flask volume 200/250 mL, used with SK330.2[13pcs] SK180.2[8pcs]	■	■
18900033	Fixing clip for flask volume 500mL, used with SK330.2[9pcs] SK180.2[8pcs]	■	■