

Portable Pressure Steam Sterilizer



Use

The product utilizes saturated steam to rapidly and efficiently sterilize articles. It can be used in hospitals, public health centers, health stations, clinics in factories and mines, scientific research institutes, etc. to sterilize medical equipment, surgical dressings, glass utensils, solutions, substrata. It can also be used as cooker in plateau areas and be used to make high quality drinking water in industrial and mining enterprises.

Characteristics

- ◆ Fully stainless steel structure.
- ◆ Hand wheel type of quick-open door structure.
- ◆ Digital display of working status, touch type key.
- ◆ Automatically shut off with beep reminding after sterilization.
- ◆ Electric heated.
- ◆ Easy to operate, safe and reliable.

LED display automation



Portable Pressure Steam Sterilizer

Technical Data

| MODEL | YX-280E18 | YX-280E24 |
|--------------------------|---------------------|---------------------|
| Technical Data | 18L φ280 × 260mm | 24L φ280 × 390mm |
| Working pressure | 0.14~0.16Mpa | |
| Working temperature | 126℃ | |
| Temperature adjust range | 105~126℃ | |
| Timer range | 0~99min | |
| Max. Safety pressure | 0.165Mpa | |
| Heat average | ≤ ± 1℃ | |
| Power | AC220V.50Hz / 2KW | |
| Transport dimension | 420 × 420 × 630 | 420 × 420 × 750 |
| G.W/N.W | 14 / 12kg | 15 / 14kg |